

## Certificate of Analysis

---

Product Name: **ATP Enumeration Reagent Buffer**

Lot Number: **390924**

Storage Conditions: **-20°C**

Expiration Date: **10/13/22**

**FOR RESEARCH USE ONLY**

**CAUTION: Not intended for human or animal therapeutic use.**

### Quality Control

---

Test	Specification	Results
Signal Stability	Signal ½ life > 5hr	Pass
Sensitivity	$RLU_{0.01\mu MATP} > RLU_{Bkg} + 2(S.D._{0.01\mu MATP} + S.D._{Bkg})$	Pass
Batch to Batch Consistency	Comparison to a Dynex reference plate	Pass

---

Manufactured and tested for Preferred Cell Systems™ by Supplier

---

Paul Bork, QA Manager

Date: